

09

TRANSFER PIPETTE

© Transfer Pipette

REF: 141

71-74



TEL:

+972-72-2786574



ADDRESS:

Bareket St. no. 7,
Caesarea North Industrial Park, Israel



E-MAIL:













sales@miniplast.co.il



WEB:

www.mini-plast.com

ICONS


	Code No.		Qty in bag
	Material		Qty in box
	Volume capacity		Qty in case for glass tube
	Additive		Qty in case for PET tube
	Drop volume		Qty in case
	Outer dimension		Sterile by irradiation
	Colour		Sterile by ethylene oxide
	Specification		Detailed information
	Temperature range		

MATERIALS

PE =	Polyethylene	COC =	Cycloolefin copolymers
PET =	Polyethylene terephthalate	HIPS =	High impact polystyrene
HDPE =	High Density Polyethylene	TPE =	Thermoplastic Elastomer
LDPE =	Low Density Polyethylene	GPPS =	General purpose polystyrene
PP =	Polypropylene	PC =	Polycarbonate
PS =	Polystyrene	P.O.M. =	Polyoxymethylene
PVC =	Polyvinyl chloride		
ABS =	Acrylonitrile Butadiene Styrene		


Transfer Pipette

- Made of superior FDA certified LDPE material, Non-toxic, adapted to draw, transfer or carry for small volume liquid.
- Optimizing process on surface tension, easy for the liquid flowing.
- High transparency, easy to observe.
- Can be bend with a certain angle, which is convenient for drawing or adding liquid in irregular or micro container.
- High elasticity, adapted to fast liquid transfer without leakage.
- Convenient and accurate for use with good repeatability.
- Heat-seal on the tip of the pipette can achieve liquid carry.
- Available in bulk or individual pack.
- Available in E.O. or Gamma radiation sterile.



REF	Original REF.	MAT						
141-001	KJ601	LDPE	65mm	10 μ L	2000		20000	Fine tip STERILE EO or STERILE R
141-002	KJ601-1	LDPE	70mm	15 μ L	2000		20000	Fine tip STERILE EO or STERILE R
141-003	KJ603	LDPE	80mm	20 μ L	2000		20000	Fine tip STERILE EO or STERILE R
141-004	KJ603-1	LDPE	103mm	12 μ L	200		12000	Fine tip STERILE EO or STERILE R
141-005	KJ604	LDPE	148mm	15 μ L	500		5000	Fine tip STERILE EO or STERILE R
141-006	KJ602	LDPE	67mm	20 μ L	2000		20000	STERILE EO or STERILE R
141-007	KJ605	LDPE	65mm	48 μ L	2000		20000	STERILE EO or STERILE R
141-008	KJ606	LDPE	67mm	48 μ L	2000		20000	STERILE EO or STERILE R
141-009	KJ607	LDPE	85mm	40 μ L	2000		20000	STERILE EO or STERILE R
141-010	KJ608	LDPE	115mm	60 μ L	2000		20000	STERILE EO or STERILE R

Transfer Pipette



141-011
(KJ609)

141-012
(KJ612)

141-013
(KJ633)

141-014
(KJ617)

141-015
(KJ618)

141-016
(KJ619)

141-017
(KJ620)


141-018
(KJ621)


141-019
(KJ622)

141-020
(KJ626)

REF	Original REF.	MAT						
141-011	KJ609	LDPE	145mm	40µL	400		4000	STERILE EO OF STERILE R
141-012	KJ612	LDPE	130mm	40µL	2000		20000	STERILE EO OF STERILE R
141-013	KJ633	LDPE	160mm	40µL		500	5000	STERILE EO OF STERILE R
141-014	KJ617	LDPE	143mm	42µL		500	5000	Graduated to 1mL STERILE EO OF STERILE R
141-015	KJ618	LDPE	145mm	45µL		500	5000	Graduated to 1mL STERILE EO OF STERILE R
141-016	KJ619	LDPE	160mm	42µL		500	5000	Graduated to 1mL STERILE EO OF STERILE R
141-017	KJ620	LDPE	150mm	48µL		500	5000	Graduated to 2mL STERILE EO OF STERILE R
141-018	KJ621	LDPE	160mm	48µL		500	5000	Graduated to 3mL STERILE EO OF STERILE R
141-019	KJ622	LDPE	162mm	52µL		500	6000	Graduated to 3mL STERILE EO OF STERILE R
141-020	KJ626	LDPE	178mm	48µL		500	5000	Graduated to 3mL STERILE EO OF STERILE R

Transfer Pipette



REF	Original REF.	MAT						
141-021	KJ623	LDPE	64mm	25µL	2000		20000	Double-bulb <input type="checkbox"/> STERILE EQ or <input type="checkbox"/> STERILE R
141-022	KJ624	LDPE	85mm	30µL	2000		20000	Double-bulb <input type="checkbox"/> STERILE EQ or <input type="checkbox"/> STERILE R
141-023	KJ645-1	LDPE	110mm	35µL	2000		20000	Double-bulb <input type="checkbox"/> STERILE EQ or <input type="checkbox"/> STERILE R
141-024	KJ641	LDPE	138mm	30µL		500	9000	<input type="checkbox"/> STERILE EQ or <input type="checkbox"/> STERILE R
141-025	KJ640	LDPE	135mm	30µL	200		12000	<input type="checkbox"/> STERILE EQ or <input type="checkbox"/> STERILE R
141-026	KJ629	LDPE	70mm	5µL	2000		20000	<input type="checkbox"/> STERILE EQ or <input type="checkbox"/> STERILE R
141-027	KJ630	LDPE	75mm	50µL	1000		10000	<input type="checkbox"/> STERILE EQ or <input type="checkbox"/> STERILE R
141-028	KJ625	LDPE	78mm	35µL	2000		20000	<input type="checkbox"/> STERILE EQ or <input type="checkbox"/> STERILE R
141-029	KJ625-1	LDPE	110mm	40µL	2000		20000	<input type="checkbox"/> STERILE EQ or <input type="checkbox"/> STERILE R
141-030	KJ606-1	LDPE	80mm	25µL		2000	20000	<input type="checkbox"/> STERILE EQ or <input type="checkbox"/> STERILE R
141-031	KJ646	LDPE	101mm	35µL		2000	20000	<input type="checkbox"/> STERILE EQ or <input type="checkbox"/> STERILE R

Transfer Pipette



REF	Original REF.	MAT						
141-032	KJ643	LDPE	155mm	40µL	500	-	5000	STERILE EQ or STERILE R
141-033	KJ644	LDPE	140mm	35µL	500	-	5000	STERILE EQ or STERILE R
141-034	KJ642	LDPE	210mm	30µL	250	-	2000	Flexible stem for tying and seal STERILE EQ or STERILE R
141-035	KJ638	LDPE	230mm	40µL	-	-	5000	STERILE EQ or STERILE R
141-036	KJ616	LDPE	350mm	55µL	200	500	2000	STERILE EQ or STERILE R
141-037	KJ620-1	LDPE	155mm	35µL(2mL)	-	-	5000	STERILE EQ or STERILE R
141-038	KJ647	LDPE	155mm	50µL	-	500	5000	STERILE EQ or STERILE R
141-039	-	LDPE	205mm	40µL(5mL)	-	500	4000	STERILE EQ or STERILE R
141-040	-	LDPE, Driper: Glass	73.4mm	20µL	2000	400	20000	STERILE EQ or STERILE R
141-041	-	LDPE, Driper: Glass	70.4mm	10µL	2000	-	20000	STERILE EQ or STERILE R