

10

PIPETTE TIP

© Pipette Tip	REF: 121	76-79
© Pipette Tip Box	REF: 122	80
© Pipettor & Reagent Reservoir	REF: 125 / 126 / 127	81
© Universal Pipette Tip	REF: 123	82-83
© Pipette tip for Automatic System	REF: 124	84-86
Reloaded tips for TECAN		84
Reloaded tips for ADDCARE		84
Boxed tips for TECAN		85
Boxed tips for AGILENT		85
Boxed tips for BECKMAN		85
Boxed tips for HAMILTON		86



TEL:

+972-72-2786574



ADDRESS:

Bareket St. no. 7,
Caesarea North Industrial Park, Israel



E-MAIL:


















sales@miniplast.co.il



WEB:

www.mini-plast.com

ICONS


	Code No.		Qty in bag
	Material		Qty in box
	Volume capacity		Qty in case for glass tube
	Additive		Qty in case for PET tube
	Drop volume		Qty in case
	Outer dimension		Sterile by irradiation
	Colour		Sterile by ethylene oxide
	Specification		Detailed information
	Temperature range		

MATERIALS

PE =	Polyethylene	COC =	Cycloolefin copolymers
PET =	Polyethylene terephthalate	HIPS =	High impact polystyrene
HDPE =	High Density Polyethylene	TPE =	Thermoplastic Elastomer
LDPE =	Low Density Polyethylene	GPPS =	General purpose polystyrene
PP =	Polypropylene	PC =	Polycarbonate
PS =	Polystyrene	P.O.M. =	Polyoxymethylene
PVC =	Polyvinyl chloride		
ABS =	Acrylonitrile Butadiene Styrene		


Pipette Tip


- ◆ The pipette tip is made of medical grade polypropylene (PP) raw material. Excellent molding process guarantee each tips straight, no bending and deformation, with high quality and precision.
- ◆ The inner surface of the low residue tip is treated by coating process, which has super hydrophobicity. Almost no sample liquid retains in the pipette tip during liquid release.
- ◆ Filtered tips prevent aerosol cross-contamination between pipette and sample.
- ◆ The optimized shape of the medical tips is designed for versatility and is perfectly adapted to pipettes like Gilson, TECAN, Eppendorf, Thermo (Finn), Kveen, Hamilton, DLab and other well-known brands at home and abroad.
- ◆ General grade and PCR grade are available according to different requirements. Produced in 100,000 class clean room, up to PCR clean grade, Free of DNAase, RNAase, pyrogen.
- ◆ Bags and boxes both available for package.



REF	Original REF.	MAT							
121-01801	KJ401	PP	-	1000μL	500	-	10000	Tip for Eppendorf	
121-01601	KJ402	PP	-	300μL	1000	-	20000	Tip for Eppendorf	
121-01401	KJ403	PP	-	200μL	1000	-	20000	Tip for Eppendorf	
121-01101	KJ404	PP	-	10μL	1000	-	50000	Tip for Eppendorf	
121-01802	KJ417-1	PP	-	1000μL	-	100	100×50	Boxed filter tips for Eppendorf	STERILE EO Or STERILE R
121-01402	KJ403-1	PP	-	200μL	-	96	96×50	Boxed filter tips for Eppendorf	STERILE EO Or STERILE R
121-01102	KJ418-1	PP	-	10μL	-	96	96×50	Boxed filter tips for Eppendorf	STERILE EO Or STERILE R






Pipette Tip




REF	Original REF.	MAT						
121-02801 (KJ405)	KJ405	PP	Ø8×71mm	1000µL	500	-	10000	Universal tip
121-02401 (KJ406)	KJ406	PP	Ø6×52mm	200µL	1000	-	20000	Universal tip
121-02101 (KJ408)	KJ408	PP	Ø5×32mm	10µL	1000	-	50000	Universal tip
121-02802 (KJ417)	KJ417	PP	Ø8×71mm	1000µL	-	100	100×50	Boxed universal filter tips STERILE EO Or STERILE R
121-02402 (KJ418)	KJ418	PP	Ø6×52mm	200µL	-	96	96×50	Boxed universal filter tips STERILE EO Or STERILE R
121-02102 (KJ418-2)	KJ418-2	PP	Ø5×32mm	10µL	-	96	96×50	Boxed universal filter tips STERILE EO Or STERILE R
121-03801 (KJ412)	KJ412	PP	Ø8×71mm	1000µL	500	-	10000	Tip for Finn
121-03401 (KJ413)	KJ413	PP	Ø6×52mm	200µL	1000	-	20000	Tip for Finn, with graduation
121-03402 (KJ414)	KJ414	PP	Ø6×52mm	200µL	1000	-	20000	Tip for Finn
121-03101 (KJ415)	KJ415	PP	Ø5×32mm	10µL	1000	-	50000	Tip for Finn
121-03802 (KJ412-1)	KJ412-1	PP	Ø8×71mm	1000µL	-	100	100×50	Boxed filter tips for Finn STERILE EO Or STERILE R
121-03403 (KJ413-1)	KJ413-1	PP	Ø6×52mm	200µL	-	96	96×50	Boxed filter tips for Finn, with graduation STERILE EO Or STERILE R
121-03404 (KJ414-1)	KJ414-1	PP	Ø6×52mm	200µL	-	96	96×50	Boxed filter tips for Finn STERILE EO Or STERILE R
121-03102 (KJ415-1)	KJ415-1	PP	Ø5×32mm	10µL	-	96	96×50	Boxed filter tips for Finn STERILE EO Or STERILE R





Pipette Tip



REF	Original REF.	MAT					
121-04501	KJ409	PP	Ø4×49mm	250µL	1000	20000	Tip for Quijing
121-04502	KJ410	PP	Ø6×50mm	250µL	1000	20000	Tip for Quijing
121-04801	KJ411	PP	Ø7.4×70mm	1000µL	500	10000	Tip for Quijing
121-14601	KJ423	PP		300µL	1000	20000	Tip for Kodak Multi-use
121-15401	KJ424	PP		200µL	1000	20000	Tip for Finn (Extended)
121-05901	KJ420	PP		5mL	300	3000	Tip for D-LAB
121-03901	KJ421	PP		5mL	300	3000	Tip for Finn
121-06901	KJ422	PP		5mL	300	3000	Tip for Gilson / Keeven
121-15901	KJ437	PP		5mL	300	3000	Multi-use
121-15902	KJ438	PP		5mL	300	3000	Multi-use



Pipette tip



REF	Original REF.	MAT				
121-12801	KJ416	PP	1000µL	500	10000	Tip for MAL
121-13801	KJ416-3	PP	1000µL	500	10000	Tip for OXFORD
121-12401	KJ416-1	PP	200µL	1000	20000	Tip for MAL
121-15402	KJ416-2	PP	200µL	1000	20000	Multi-use
121-05101	KJ430	PP	10µL	1000	50000	Tip for D-LAB
121-05401	KJ419-3	PP	200µL	1000	20000	Tip for D-LAB
121-05081	KJ430-2	PP	1000µL	500	10000	Tip for D-LAB
121-15403	KJ427	PP	200µL	1000	20000	Multi-use
121-15404	KJ428	PP	200µL	1000	20000	Multi-use
121-15405	KJ428-1	PP	200µL	1000	20000	Multi-use
121-15406	KJ428-2	PP	200µL	1000	20000	Multi-use
121-15801	KJ429	PP	1000µL	500	10000	Multi-use
121-15802	KJ429-1	PP	1000µL	500	10000	Multi-use
121-15803	KJ429-2	PP	1000µL	500	10000	Multi-use
121-15804	KJ429-3	PP	1000µL	500	10000	Multi-use
121-15301	KJ434	PP	100µL	1000	20000	Multi-use

Pipette Tip Box



REF	Original REF.	MAT		
122-16401	KJ447	PP	200µL	50Box
122-16101	KJ448	PP	10µL	50Box
122-16801	KJ411-4	PP	1000µL	50Box
122-16901	KJ422-1	PP	5mL	30Box
122-16001	-	PP	10mL	30Box




Pipettor & Reagent Reservoir


 125-0001
 (KJ411-1)

 125-0002
 (KJ449)


126-0200

 127-0001
 (KJ426)

REF	Original REF.	MAT			
125-0001	KJ411-1	PVC	250mL	800	Reagent reservoir
125-0002	KJ449	PP	10mL×10	Customized	10-Channel Reagent reservoir
126-0200	-	ABS		Customized	Thermofisher / BIO-DL Pipettor, 1-channel / 8-channel / 12-channel
127-0001	KJ426	-		Customized	Stand for Pipettor

Universal Pipette Tip [Pyrogen-free]

- ◆ The pipette tips are made of medical-grade polypropylene by advanced production techniques. Straight without any bending or deformation, ensuring high quality and handling precision.
- ◆ Produced in a ISO Class 6 cleanroom, these pipette tips meet high cleanliness requirements to ensure superior purity and quality, and are free of DNase, RNase, and pyrogens.
- ◆ **MINIPLAST** pipette tips feature a universally compatible design, perfectly fitting pipettors from a variety of well-known brands.
- ◆ The tips with filters prevent aerosol cross-contamination between the pipettor and samples.
- ◆ Available in multiple packaging options, including bulk pack bags or specialized tip boxes.

10uL Universal Tip

10uL Universal Tip, Extended



200uL Universal Tip

200uL Universal Tip, Extended










1000uL Universal Tip

1000uL Universal Tip, Extended



Universal Pipette Tip

REF							
123-15101	10µL	32mm	PP	1000		50000	Clear, Bagged
123-15102	10µL	32mm	PP		96	4800	Clear, Boxed
123-15103	10µL	32mm	PP	1000		50000	Clear with Filter, Bagged
123-15104	10µL	32mm	PP		96	4800	Clear with Filter, Boxed
123-15105	10µL (Extended)	45mm	PP	1000		50000	Clear, Bagged
123-15106	10µL (Extended)	45mm	PP		96	4800	Clear, Boxed
123-15107	10µL (Extended)	45mm	PP	1000		50000	Clear with Filter, Bagged
123-15108	10µL (Extended)	45mm	PP		96	4800	Clear with Filter, Boxed
123-15401	200µL	51mm	PP	1000		20000	Yellow or Clear, Bagged
123-15402	200µL	51mm	PP		96	4800	Yellow or Clear, Boxed
123-15403	200µL	51mm	PP	1000		20000	Yellow or Clear with Filter, Bagged
123-15404	200µL	51mm	PP		96	4800	Yellow or Clear with Filter, Boxed
123-15405	200µL (Extended)	60mm	PP	1000		20000	Clear, Bagged
123-15406	200µL (Extended)	60mm	PP		96	4800	Clear, Boxed
123-15407	200µL (Extended)	60mm	PP	1000		20000	Clear with Filter, Bagged
123-15408	200µL (Extended)	60mm	PP		96	4800	Clear with Filter, Boxed
123-15801	1000µL	71mm	PP	500		10000	Clear or Blue, Bagged
123-15802	1000µL	71mm	PP		100	5000	Clear or Blue, Boxed
123-15803	1000µL	71mm	PP	500		10000	Clear or Blue with Filter, Bagged
123-15804	1000µL	71mm	PP		100	5000	Clear or Blue with Filter, Boxed
123-15805	1000µL (Extended)	102mm	PP	500		10000	Clear, Bagged
123-15806	1000µL (Extended)	102mm	PP		100	5000	Clear, Boxed
123-15807	1000µL (Extended)	102mm	PP	500		10000	Clear with Filter, Bagged
123-15808	1000µL (Extended)	102mm	PP		100	5000	Clear with Filter, Boxed

Pipette tip for Automatic System

- These tips are designed for excellent liquid treatment result, Under strict quality control procedures to ensure consistency to get excellent tip straightness and reduced variation coefficient.
- The tips are produced in Class 100,000 clean workshop, made from imported polypropylene material. Each batch is strictly evaluated, DNase-free, RNase-free, pyrogen-free, Without endotoxin or PCR inhibitor.

Reloaded tips for TECAN



Reloaded tips for ADDCARE



REF			Paper box		Inner box		Outer carton	
124-07401	200µL	127×88×95mm	96+96 PCS	[96+96]×12	[96+96]×12×4	For TECAN automatic system 96 PCS Clear tip + 96 PCS Conductive Tip		
124-07801	1000µL	127×88×159mm	96+96 PCS	[96+96]×12	[96+96]×12×4	For TECAN automatic system 96 PCS Clear tip + 96 PCS Conductive Tip		
124-10601	300µL	125×82×96mm	192 PCS	192×12	192×12×4	For ADDCARE automatic system 96+96 PCS Clear tip		
124-10701	800µL	125×82×126mm	192 PCS	192×12	192×12×4	For ADDCARE automatic system 96+96 PCS Clear tip		

Pipette tip for Automatic System

Boxed tips for TECAN


 124-07402
 (KJ441)

 124-07802
 (KJ442)

 124-07403
 (KJ443)





 124-07803
 (KJ444)

Boxed tips for AGILENT


 124-08201
 (KJ445)

Boxed tips for BECKMAN






 124-09501
 (KJ446)

REF	Original REF.				
124-07402	KJ441	200µL	96	50	Boxed tips for TECAN
124-07802	KJ442	1000µL	96	50	Boxed tips for TECAN
124-07403	KJ443	200µL	96	50	Boxed tips for TECAN
124-07803	KJ444	1000µL	96	50	Boxed tips for TECAN
124-08201	KJ445	30µL	384	50	Boxed tips for Agilent
124-09501	KJ446	250µL	96	50	Boxed tips for Beckman

Pipette tip for Automatic System

Boxed tips for HAMILTON



REF				
124-17051	50µL	96	50	50ul Hamilton conductive tip (Customized clear tip)
124-17601	300µL	96	50	300ul Hamilton conductive tip (Customized clear tip)
124-17602	300µL	96	50	300ul Hamilton conductive tip, extended (Customized clear tip)
124-17801	1000µL	96	50	1000ul Hamilton conductive tip
124-17802	1000µL	96	50	1000ul Hamilton clear tip